



Shaver Blades Einsetzbar für Shaver Systeme:

Shaver Blades Fitting for Shaver Systems:



Tiefenmarkierung 50mm depth marking

REMA "exact-line"

		autoclavable				sterile single use		
	Farbe Color	Ø 2,5 mm L: 80 mm	Ø 3,4 mm L: 130 mm	Ø 4,2 mm L: 130 mm	Ø 5,5 mm L: 130 mm	Ø 3,4 mm L: 130 mm	Ø 4,2 mm L: 130 mm	Ø 5,5 mm L: 130 mm
Aggressive Full-Radius Resector 		...210	...310	...410	...510	...310S	...410S	...510S
Double Aggressive Full-Radius 				...415	...515			
Meniscus cutter 			...320	...420	...520	...320S	...420S	...520S
End Cutter 				...425	...525			
Full-Radius Resector 		...230	...330	...430	...530	...330S	...430S	...530S
Incisor 			...340	...440	...540	...340S	...440S	...540S
Cutter 				...450				
New Round Burr 				...460	...560			...560S
Tapered Burr 					...570			...570S
New Oval Burr 				...480	...580			...580S

Verpackungseinheit = 5 Stück Packaging unit = 5 pieces



Cleaning-Adapter

REINIGUNGSBÜRSTE PACK 10 STCK 24-00016

CLEANING BRUSH PACKED 10 EACH 24-00016

FDA K060172

