



REMA[®]

SURGERY



Made in Germany

5 THINGS IMPORTANT about **Shaver Blades** by REMA

02

The REMA rotatory blades and burs are built of **stainless steel**, the couplings are made of **anodized aluminium** or **plastic PPSU medical grade**.

04

Available in reusable and steril versions,
AUTOCCLAVABLES 134°C.

01

With **8 different couplings** compatibles with another shaver brands.



03

Wide range of geometric forms.



05

There are different lengths, straight and curves, **designed for small joints**, knee, shoulder and hip.

Small joints

80 mm Straight
ø 2.5 mm



Knee and shoulder

130 mm Straight
ø 3.4 mm, 4.2 mm y 5.5 mm



130 mm Curve
ø 4.2 mm y 5.5 mm



Hip

210 mm Straight
ø 3.4 mm, 4.2 mm y 5.5 mm

